DOSYVVEL®

PERISTALTIC PUMP

Small and efficient





PERISTALTIC PUMP

Small and efficient

DOSYWEL®

Quick and easy calibration

Calibration reminder at the beginning of the cycle or during the distribution Up to 1% of precision/accuracy





Compact and light

Easily transportable with its handle Small footprint

3.3 kg



3 ways to distribute

Fingerstop design

Avoid cross-contamination





| 9ml | | 1 | second |
|-------|-------------------|----|--------|
| 18ml | - | 3 | second |
| 225ml | \longrightarrow | 10 | second |



IT WILL FOLLOW YOU ANYWHERE!













TECHNICAL DATA SHEET Manual, automatic or with footswitch Distribution 0.1 to 9999 mL Adjustable volume Adjustable speed 1-130% - 10-550 rpm Easy and fast at the beginning of a cycle Calibration and/or each X doses Pause time 0.1 to 99 seconds Left or right Rotation way Programs USB port and Ethernet Connectivity 1.6mm / 3.2mm / 4.8mm / 6.4mm Internal tube dimensions $\mathsf{Versilion}^{\mathsf{TM}} - \mathsf{rotor} \ \mathsf{Pumpsil}^{\ \mathsf{TM}}$ Dimensions (W x H x D) 250 x 190 x 210 mm Weight 3.3 kg Power supply 100-250v - 50/60 Hz - 90W Warranty 1 year Country of origin France - European Union

Choose the good set of tubing according to your needs:

| Tube Ø (mm) | 9ml | 18ml | 225ml | Flow(ml/min) | Repeatability (± ml) |
|------------------|------------|------------|------------|--------------|----------------------|
| Ø 1.6 (DOSW2006) | ± 0.037 ml | ± 0.037 ml | ± 0.037 ml | 102 | 0,007 |
| Ø 3.2 (DOSW2003) | ± 0.159 ml | ± 0.159 ml | ± 0.159 ml | 406 | 0.055 |
| Ø 4.8 (DOSW2002) | ± 0.365 ml | ± 0.365 ml | ± 0.365 ml | 936 | 0.06 |
| Ø 6.4 (DOSW2001) | ± 0.53 ml | ± 0.53 ml | ± 0.53 ml | 1271 | 0.07 |

SPECIFICATIONS

OPTIONS:

- Additional pump
- Footswitch
- Ticket printer
- Keyboard

FOOTSWITCH



Visit our website:

alliance-bio-expertise.com



Delivered with:

DOSYWEL (DOSW1001) - Set of tubing 3.2mm with distribution nozzle (DOSW2003) - Power supply - USB key with user manual and quality report

Accessories:

DILW2001 - Additional pump

DOSW2001 - Set of tubing 6.4mm with dispensing nozzle

DOSW2002 - Set of tubing 4.8mm with dispensing nozzle

DOSW2003 - Set of tubing 3.2mm with dispensing nozzle

DOSW2006 - Set of tubing 1.6mm with dispensing nozzle

DILW3003 - Precut rotor tubing 3.2mm - 5 pcs (available in 6.4mm, 4.8mm and 1.6mm)

DILW3009 - Bulk rotor tubing 3.2mm - 5 meters (available in 6.4mm, 4.8mm and 1.6mm)

DOSW2004 - Classical distribution nozzle

DOSW2008 - Precision distribution nozzle

DOSW2005 - Stainless steel tube holder 12 positions

DILW2006 - Priming nozzle for broth bottles

DILW2007 - Trocar (connector) for broth bags

Qualification:

DOSW4000 - DOSYWEL IQ/OQ

Connectivity:

DILW2011 - Ticket printer

DILW2013 - Footswitch

DILW2038 - Clean Keyboard - Tactile wireless keyboard

Ready to use broth bags:

BGWR2001 - Support stand for BAGGYWEL (Broth bags)

BGWR3001 - BAGGYWEL Buffered Peptone Water / broth bags $2 \times 5L$ (also available in $3 \times 3L$)

BGWR3005 - BAGGYWEL Half Fraser / broth bags 2 x 5L (also available in 3 x 3L)







A FULLY COMMITTED TEAM

OF SERVICE ENGINEERS AND MICROBIOLOGISTS





A CONNECTED EQUIPMENT GENERATION FOR TRACEABILITY AND PRODUCTIVITY








FOOD SAMPLE PREPARATION

AUTOMATED MEDIA PREPARATION

SAMPLE COLLECTION














BOOK YOUR DEMO!

https://www.alliance-bio-expertise.com/en/store-locator/











ZI de Courbouton, Le Tremplin - 35480 GUIPRY - FRANCE

Phone: +33(0)2.40.51.79.53■ alliance-bio-expertise.com